

# Disinfection of Hybrid Contact Lenses

---

This material will help you understand how to disinfect your hybrid contact lenses.

## **Why is it important to disinfect my contact lenses?**

When you wear contacts, materials such as mucous, protein, and other debris can build up on the surface during the day. Disinfecting your lenses removes these materials as well as harmful germs. This helps prevent you from developing an eye infection or another serious complication.

## **How often do I need to disinfect my contact lenses?**

You should clean your contact lenses every time you take them out of your eyes. It is best to get into the routine of cleaning them immediately after removal.

## **How do I disinfect my contact lenses?**

When you are disinfecting your hybrid lenses, it is important to follow these steps:

1. Wash your hands with soap and water. Dry your hands with a lint-free towel.
2. After removing the lens from your eye, place it on your palm so it sits like a bowl.
3. Hybrid lenses can be cleaned with a hydrogen peroxide based solution such as Clear Care® or with a multi-purpose solution for soft lenses.
4. If you are using clear Care®, place lenses into lens holders. Then, thoroughly rinse the lenses for 5 seconds through the basket. Fill the

case with fresh Clear Care® solution and tighten the cap. Allow the lenses to soak for at least 6 hours.

5. If you are using a multi-purpose solution for soft lenses, place two drops of the solution recommend by your eye doctor on the lens. Gently rub the lens back and forth with the pad of your finger. Rinse the lens with a preservative free saline solution recommended by your doctor. Fill the case with fresh multi-purpose solution and allow it to soak for at least 6 hours.

**For more information, scan these codes with your smartphone or visit the websites listed.**



<http://www.synergeyes.com/consumer/ultrahealth/handling-care/>

Disclaimer: This document contains information and/or instructional materials developed by the University of Michigan Health System (UMHS) for the typical patient with your condition. It may include links to online content that was not created by UMHS and for which UMHS does not assume responsibility. It does not replace medical advice from your health care provider because your experience may differ from that of the typical patient. Talk to your health care provider if you have any questions about this document, your condition or your treatment plan.

Author: Kristin Maurer, MPH candidate  
Reviewers: Amy Lagina, OD,FAAO and Gale Oren, MILS

Unless otherwise noted, Patient Education by [University of Michigan Health System](#) is licensed under a [Creative Commons Attribution-NonCommercial-ShareAlike 3.0 Unported License](#).

Last Revised 4/2015